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Human platelet concentrates.

Presentation:

Vial containing approximately 100 µg of PF4 in the tetrameric form.

Reconstitution:

Each vial must be restored with 1 ml distilled water. When restored, the solution contains a 100 µg/ml solution of PF4 in a 0.15 M sodium chloride isotonic buffer. The vial can be restored with 0.25 ml or 0.50 ml distilled water for keeping a sodium chloride concentration of 0.6 M or 0.3 M. The PF4 concentration is then of 400 µg/ml or 200 µg/ml.

Nota : A high salt concentration (≥ 0.3 M) is recommended for long term storage (> 6h) after reconstitution.

Excipients:

– Glycine and hepes buffer - 15 mg sodium chloride/vial

Purity:

One major band of 8,000 daltons on SDS-PAGE..

Activity:

Tested in a clotting assay for the heparin inhibition: 1 vial (100 µg) neutralizes at least 1.25 IU of heparin.

Viral safety:

The material used for PF4 preparation was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Use:

For in vitro use only.
All research studies and protocols where a source of highly purified PF4 is necessary.

Stability:

Restored : - 8 hours at room temperature
- 24 hours at 2-8 °C
- 2 months frozen at -30°C or below

Lyophilized : Until the expiration date printed on the vial.

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