

CERTIFICATE OF ANALYSIS

CHROMOGENIX

PRODUCT:	COATEST APC Resistance Control Plasma Level 1	PRODUCT # :	82265063
EXP. DATE:	2006-03	LOT # :	N0347918
VENDOR:	INSTRUMENTATION LABORATORY	ASSAYED PER QCM:	C00423

TEST DESCRIPTION	SPECIFICATIONS		UNITS	RESULTS	DISPOSITION	
	MIN	MAX				
Analysis with COATEST APC R:						
APT Time	<u>CaCl₂</u>	28	34	seconds	31	PASS
	<u>APC/CaCl₂</u>	90	--	seconds	101	PASS
APC Ratio		2.7	4.5	--	3.3	PASS

Analysis with COATEST APC R V:						
APT Time	<u>CaCl₂</u>	--	38	seconds	36	PASS
	<u>APC/CaCl₂</u>	75	--	seconds	103	PASS
APC Ratio		2.3	3.0	--	2.9	PASS

Kit Components:

DESCRIPTION	LOT NUMBER	EXP. DATE	VOLUME
Control Plasma Level 1	N0347918	2006-03	1 ml (x 5)

We hereby certify that the above product(s):

- have been manufactured and tested according to the requirements of their Manufacturing and Quality Control procedures and have been found to meet all conditions set forth with these documents;

PREPARED BY: D. Tencza

DATE: 05/24/04

APPROVAL SIGNATURE: J. Facciola

DATE: 05/24/04