

CERTIFICATE OF ANALYSIS

**CHROMOGENIX**

PRODUCT:	COATEST APC Resistance Control Plasma Level 1	PRODUCT # :	82265063
EXP. DATE:	Nov-05	LOT # :	N1136040
VENDOR:	INSTRUMENTATION LABORATORY	ASSAYED PER QCM:	C00423

TEST DESCRIPTION	SPECIFICATIONS			RESULTS	DISPOSITION
Analysis with COATEST APC R:	MIN	MAX	UNITS		PASS/FAIL
APT Time <u>CaCl<sub>2</sub></u>	28	34	seconds	30	Pass
<u>APC/CaCl<sub>2</sub></u>	90	--	seconds	92	Pass
APC Ratio	2.7	4.5	--	3.1	Pass

Analysis with COATEST APC R V:					
APT Time <u>CaCl<sub>2</sub></u>	--	38	seconds	33	Pass
<u>APC/CaCl<sub>2</sub></u>	75	--	seconds	95	Pass
APC Ratio	2.3	3.0	--	2.8	Pass

**Kit Components:**

DESCRIPTION	LOT NUMBER	EXP.DATE	VOLUME
Control Plasma Level 1	N1136040	2005-11	1 ml (x 5)

We hereby certify that the above product(s):

- have been manufactured and tested according to the requirements of their Manufacturing and Quality Control procedures and have been found to meet all conditions set forth with these documents;

PREPARED BY: M. Brunswick

DATE: 12/22/03

APPROVAL SIGNATURE: R. Bustamante

DATE: 12/22/03