

**REAADS Coagulation  
Controls 1 and 2**  
*For In Vitro Diagnostic Use*

**REAADS Coagulation  
Controls 1 and 2**

**REAADS Coagulation  
Control 1**

Product #: 701-001  
(Box of 10 vials to be reconstituted to 1ml)

**REAADS Coagulation Control 2**

Product #: 702-001  
(Box of 10 vials to be reconstituted to 1ml)

**Product Description**

REAADS Coagulation Controls are assayed control plasma for *in vitro* diagnostic use in coagulation studies. The controls are lyophilized and to be reconstituted in 1.0 mL deionized water.

REAADS Coagulation Control 1 serves as a normal control. REAADS Coagulation Control 2 serves as a abnormal control.

**Principle**

REAADS Coagulation Controls are tested in the same manner as citrated patient plasma samples to assess the performance of each assay run, for the parameters and by the methods listed on the back page.

**REAADS Coagulation  
Controls 1 and 2**  
*For In Vitro Diagnostic Use*

<b>Test</b>	<b>Procedure</b>
Protein C	REAADS ELISA
Protein S, Total	REAADS ELISA
Protein S, Free	REAADS ELISA / PEG method
Protein S, Free Monoclonal	REAADS Monoclonal ELISA
von Willebrand Factor Antigen	REAADS ELISA
Prothrombin Time	One-stage method
Partial Thromboplastin	One-stage method
Fibrinogen	Functional (Clauss)
Factor II	One-stage method
Factor V	One-stage method
Factor VII	One-stage method
Factor VIII	One-stage method
Factor VIII	Kinetic / endpoint
Factor IX	One-stage method
Factor X	One-stage method
Factor X	Kinetic / endpoint
Factor XI	One-stage method
Factor XII	One-stage method
Plasminogen	Immunological - RID
Plasminogen	Kinetic / endpoint
Antithrombin III	Immunological - RID
Antithrombin III	Kinetic / endpoint
Ristocetin Cofactor	Helena PACKS-4 aggregometer
von Willebrand Factor Antigen	Laurell Rocket EID
Protein C	Laurell Rocket EID
Protein C	Kinetic / endpoint
Protein S, Total	Laurell Rocket EID
Protein S, Free	Laurell Rocket EID / precipitation
$\alpha_2$ -Antiplasmin	Kinetic / endpoint



**Corgenix, Inc.**

12061 Tejon Street  
Westminster, Colorado 80234  
Toll Free: 800-729-5661  
Telephone: 303-457-4345  
Fax: 303-457-4519  
Email: techsupport@corgenix.com  
www.corgenixonline.com

**Corgenix (UK) Ltd.**

75 Broadway  
Peterborough, PE1 1SY UK  
Telephone: +44 (0) 1733 296800  
Fax: +44 (0) 1733 296809  
Email: sales@corgenix.co.uk  
www.corgenix.co.uk