

*Introducing Products From*

# DiaPharma

**Now Available!**  
For research use only

## Quality Research Products for the Study of Thrombosis, Hemostasis, and Autoimmunity

### DiaPharma FPA

Competitive Enzyme Immuno Assay of Fibrinopeptide A.

- Ready to use reagent complete kit
- Optimized, standardized assay
- Microplate is already precoated and stabilized
- High sensitivity and specificity
- Detection threshold < 0.5 ng/ml

FPA is a marker of hypercoagulable states induced in many pathological conditions such as DIC, SLE, malignancy, and some infections.

### DiaPharma TAFI: Ag

Enzyme Immuno Assay of TAFI (Thrombin Activatable Fibrinolysis Inhibitor) protein in human plasma or serum.

- Disposable, complete ready to use kit
- Optimized, standardized assay
- Measures the total protein
- Insensitive to plasma storage conditions
- Total assay time of only 2 – 4 hours
- Result expressed as % of normal plasma
- Assay range from 1 – 200 ng/ml



TAFI delays fibrinolysis, and high concentrations of TAFI can be a factor predisposing to thrombotic diseases.

### DiaPharma Annexin V

A new reliable assay for the measurement of Annexin V.

- Annexin V released correlates with platelet and blood cell surface lesion
- Disposable, complete ready to use kit
- Measures native or recombinant Annexin V
- No interference of rheumatoid factor
- Detection threshold 0.2 ng/ml

Annexin V is a protein present in tissue, platelets, red blood cells and leukocytes. It inhibits blood coagulation and phospholipase A, as well as pro-coagulant and pro-inflammatory activities of dying cells during apoptosis.



### KIT CONTENTS:

- 6 x 16 well microplate coated with F(ab)'2 fragments from rabbit polyclonal antibodies specific for Annexin V
- Anti-Annexin V – peroxidase immunoconjugate
- Annexin V standard
- Ready to use sample diluent
- Ready to use conjugate diluent
- Substrate tablets (H<sub>2</sub>O<sub>2</sub>/OPD)

For a complete listing of products available from DiaPharma, please call Customer Service at (800) 526-5224 or Technical Service at (800) 447-3846, or visit [www.diapharma.com](http://www.diapharma.com).

DiaPharma brand products are distributed by DiaPharma Group, Inc. in the United States and Canada.

**diaPharma**  
φ

*Introducing Products From*

# DiaPharma

## Quality Research Kits and Reagents for the Study of Thrombosis, Hemostasis, and Autoimmunity

**FOR RESEARCH USE ONLY**

- **Extensive line of ELISA kits for the exploration of hemostasis factors, fibrinolysis, antiphospholipid antibodies, blood cell and platelet markers.**
  - Ready to use and complete Enzyme Immuno assays featuring disposable 6x16 well strips
  - Easy to perform, fully optimized and reliable
  - High reproducibility and sensitivity
  
- **High quality purified human and bovine proteins and enzymes involved in coagulation and fibrinolysis**
  - Standardized, optimized, and highly reproducible
  - High purity, high specific activity, and excellent stability
  - Useful for research protocols, quality control procedures, and developing specific functional assays
  
- **Purified antibodies and conjugates**
  - High activity and avidity, excellent specificity, homogeneous reactivity to the specific protein
  - Low background noise
  - Suitable for gel precipitation techniques, immunofixation and Western blot protocols

*Distributed by:*

DiaPharma Group, Inc.	1-800-526-5224 (to order)
8948 Beckett Road	1-800-447-3846 (tech support)
West Chester, OH 45069	info@diapharma.com (e-mail)
USA	www.diapharma.com (web)

