



412 South Lafayette Blvd.
South Bend, IN 48001

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Russells' Viper Venom

Concentration: 1.00 mg/ml (before lyophilization)

Allquots: 1 x 1.00 mg

Volume/aliquot: 1.000 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ 1mM
Benzamidine/ pH 7.4

Storage: 4 °C

Expiration Date: 03/2008

PO: N/A

The Russells' Viper Venom is homogeneous as judged by SDS-PAGE and has a high specific clotting activity.

We suggest hydrating the protein to the original volume with water.

Thank You,
Quality Assurance

Sample