



Product Data Sheet for:
Sheep anti-Human Vitronectin
Affinity Purified IgG

Catalogue #: SAVN-AP

Lot #: XXXX

Expiry: XXXX

Rev : 702

Immunogen: Human vitronectin (Vn) purified from plasma. This protein is also known as S-Protein.

Presentation: 0.146 ml of immunoglobulin from sheep antiserum affinity purified on immobilized vitronectin. The IgG concentration is 3.43 mg/ml in 10 mM HEPES , pH 7.4, containing 0.15 M NaCl and 50% glycerol (v/v).

Storage: Store in the freezer between -10 to -20°C. Gently mix before use.

Specificity: This antibody is specific for vitronectin as demonstrated by immunoelectrophoresis. A single precipitin arc is observed with normal plasma. No arc was observed with vitronectin depleted plasma.

Application: This product can be used for immunostaining and in immunoassays for human vitronectin and its complexes.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

For Research Use Only

Related products:

Product: SAVN-IG	Sheep anti-human Vitronectin, purified whole IgG from antiserum	10 mg vial
Product: SAVN-AP	Sheep anti-human Vitronectin, affinity purified IgG	0.5 mg vial